

Catalog #: 52848

Aliases

TRIM45; Tripartite Motif Containing 45; RNF99; E3 Ubiquitin-Protein Ligase TRIM45; FLJ13181; Tripartite Motif-Containing 45; Ring Finger Protein 99; RING Finger Protein 99; EC 2.3.2.27

Background

Gene Name: TRIM45

NCBI Gene Entry: [80263](#)

UniProt Entry: [Q9H8W5](#)

Application Information

Molecular Weight: Predicted, 64 kDa; observed, 64 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human TRIM45

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

