

Catalog #: 52853

Aliases

PDGFB; Platelet Derived Growth Factor Subunit B; Becaplermin; SSV; SIS; Platelet-Derived Growth Factor Beta Polypeptide (Simian Sarcoma Viral (V-Sis) Oncogene Homolog); Platelet-Derived Growth Factor Beta Polypeptide; Platelet-Derived Growth Factor Subunit B; Platelet-Derived Growth Factor B Chain; Proto-Oncogene C-Sis; PDGF Subunit B; PDGF-2; PDGF2; Platelet-Derived Growth Factor, Beta Polypeptide (Oncogene SIS); Epididymis Secretory Sperm Binding Protein; Platelet-Derived Growth Factor; PDGF, B Chain; Oncogene SIS; IBGC5; C-Sis

Background

Gene Name: PDGFB

NCBI Gene Entry: [5155](#)

UniProt Entry: [P01127](#)

Application Information

Molecular Weight: Predicted, 27 kDa; observed, 27 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human PDGFB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:300

Immunocytochemistry (IC): 1:50-1:200

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data

