

Catalog #: 52855

Aliases

MYH14; Myosin Heavy Chain 14; KIAA2034; MHC16; MYH17; Myosin, Heavy Chain 14, Non-Muscle; Myosin Heavy Chain, Non-Muscle IIc; Non-Muscle Myosin Heavy Chain IIc; Myosin, Heavy Polypeptide 14; NMHC II-C; Myosin-14; FLJ13881; DFNA4; Nonmuscle Myosin Heavy Chain II-C; MYH14 Variant Protein; NMHC-II-C; FP17425; DFNA4A; Myosin; PNMHH

Background

NCBI Gene Entry: [79784](#)

UniProt Entry: [Q7Z406](#)

Application Information

Molecular Weight: Predicted, 228 kDa, observed, 228 kDa

Clonality: Rabbit polyclonal antibody

Clone ID: N/A

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunofluorescence (IF)

Immunogen

A synthesized peptide derived from human MYH14

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM