

Catalog #: 52888

Aliases

FAR1; Fatty Acyl-CoA Reductase 1; SDR10E1; MLSTD2; Short Chain Dehydrogenase/Reductase Family 10E, Member 1; Male Sterility Domain-Containing Protein; FLJ22728; Male Sterility Domain Containing 2; EC 1.2.1.84; EC 1.2.1.N2; CSPSD; PFCRD

Background

Gene Name: FAR1

NCBI Gene Entry: [84188](#)

UniProt Entry: [Q8WVX9](#)

Application Information

Molecular Weight: Predicted, 59 kDa; observed, 59 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FACR1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from Hela cells, primary antibody was diluted at 1:1,000, 4°C over night.