

Catalog #: 52893

Aliases

MAFB; MAF BZIP Transcription Factor B; KRML; V-Maf Avian Musculoaponeurotic Fibrosarcoma Oncogene Homolog B; Transcription Factor MafB; MAFB/Kreisler Basic Region/Leucine Zipper Transcription Factor; V-Maf Musculoaponeurotic Fibrosarcoma Oncogene Homolog B; Kreisler (Mouse) Maf-Related Leucine Zipper Homolog; Kreisler Maf-Related Leucine Zipper Homolog; DURS3; Maf-B; MCTO

Background

Gene Name: MAFB

NCBI Gene Entry: [9935](#)

UniProt Entry: [Q9Y5Q3](#)

Application Information

Molecular Weight: Predicted, 36 kDa; observed, 35 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human MAFB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from U2OS cells, primary antibody was diluted at 1:1,000, 4°C over night.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM