

Catalog #: 52906

Aliases

TMEM30A; Transmembrane Protein 30A; CDC50A; C6orf67; P4-ATPase Flippase Complex Beta Subunit TMEM30A; Cell Cycle Control Protein 50A; FLJ10856; Chromosome 6 Open Reading Frame 67

Background

Gene Name: TMEM30A

NCBI Gene Entry: [55754](#)

UniProt Entry: [Q9NV96](#)

Application Information

Molecular Weight: Predicted, 41 kDa; observed, 39 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TMEM30A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from DU145 cells, primary antibody was diluted at 1:1,000, 4°C over night.