

Anti-Pan Methyl Lysine Rabbit Polyclonal Antibody



Catalog #: 52908

Application Information

Clonality: Rabbit polyclonal antibody

Species Reactivity: Species independent

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of Pan Methyl Lysine

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

