

Catalog #: 52910

Aliases

NEUROG3; Neurogenin 3; BHLHa7; Math4B; Atoh5; Ngn3; Class A Basic Helix-Loop-Helix Protein 7; Protein Atonal Homolog 5; Neurogenin-3; NGN-3; BHLHA7; ATOH5; NGN3

Background

Gene Name: NEUROG3

NCBI Gene Entry: [50674](#)

UniProt Entry: [Q9Y4Z2](#)

Application Information

Molecular Weight: Predicted, 23 kDa; observed, 25 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunohistochemistry-Paraffin (IHC-P), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Neurogenin 3

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Immunohistochemistry-Paraffin (IHC-P): 1:100-1:500

Immunocytochemistry (IC): 1:10-1:50

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-Neurogenin 3 antibody (Cat#52910). Mouse liver tissue lysate were loaded and separated by SDS-PAGE. The blot was incubated with anti-Neurogenin 3 antibody (Cat#52910, 1:1,000) and secondary antibody respectively.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM