

Catalog #: 52911

Aliases

TAAR1; Trace Amine Associated Receptor 1; TAR1; TA1; Trace Amine Receptor 1; TRAR1; Trace Amine-Associated Receptor 1; TaR-1

Background

Gene Name: TAAR1

NCBI Gene Entry: [134864](#)

UniProt Entry: [Q96RJ0](#)

Application Information

Molecular Weight: Predicted, 39 kDa; observed, 37 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human TAAR1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

