

Catalog #: 52914

Aliases

FTSJ1; FtsJ RNA 2'-O-Methyltransferase 1; CDLIV; TRMT7; JM23; SPB1; TRNA (Cytidine(32)/Guanosine(34)-2'-O)-Methyltransferase; 2'-O-Ribose RNA Methyltransferase TRM7 Homolog; FtsJ RNA Methyltransferase Homolog 1; MRX44; TRM7; MRX9; Putative TRNA (Cytidine(32)/Guanosine(34)-2'-O)-Methyltransferase; TRNA Methyltransferase 7 Homolog (S. Cerevisiae); FtsJ RNA Methyltransferase Homolog 1 (E. Coli); Putative Ribosomal RNA Methyltransferase 1; RRNA (Uridine-2'-O)-Methyltransferase; TRNA Methyltransferase 7 Homolog; Mental Retardation, X-Linked 44; Mental Retardation, X-Linked 9; Protein FtsJ Homolog 1; Cell Division Protein; FtsJ-Like Protein 1; EC 2.1.1.205; XLID9

Background

Gene Name: FTSJ1

NCBI Gene Entry: [24140](#)

UniProt Entry: [Q9UET6](#)

Application Information

Molecular Weight: Predicted, 36 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FTSJ1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data

