

Catalog #: 52916

Aliases

IL20RB; Interleukin 20 Receptor Subunit Beta; IL-20R2; DIRS1; Fibronectin Type III Domain Containing 6; FNDC6; Interleukin 20 Receptor Beta Subunit; Interleukin-20 Receptor Subunit Beta; IL-20 Receptor Subunit Beta; IL-20R-Beta; MGC34923; IL-20RB; Interleukin 20 Receptor Beta; Interleukin-20 Receptor II

Background

Gene Name: IL20RB

NCBI Gene Entry: [53833](#)

UniProt Entry: [Q6UXL0](#)

Application Information

Molecular Weight: Predicted, 35 kDa; observed, 35 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human IL20RB

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions


Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:100-1:300

Immunocytochemistry (IC): 1:50-1:200

Note: This product is for research use only.

Validation Data



Western blotting analysis of various cells using IL20RB antibody.