

Catalog #: 52917

Aliases

HSPG2; Heparan Sulfate Proteoglycan 2; Perlecan; PRCAN; Basement Membrane-Specific Heparan Sulfate Proteoglycan Core Protein; Perlecan Proteoglycan; HSPG; SJS1; PLC; Schwartz-Jampel Syndrome 1 (Chondrodystrophic Myotonia); Endorepellin (Domain V Region); Endorepellin; SJA; SJS

Background

Gene Name: HSPG2
NCBI Gene Entry: [3339](#)
UniProt Entry: [P98160](#)

Application Information

Molecular Weight: Predicted, 469 kDa
Clonality: Rabbit polyclonal antibody
Species Reactivity: Human, mouse
Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human HSPG2

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

