

Catalog #: 52918

Aliases

GREM1; Gremlin 1, DAN Family BMP Antagonist; DAND2; DRM; Cysteine Knot Superfamily 1, BMP Antagonist 1; CKTSF1B1; Gremlin; HMPS; Down-Regulated In Mos-Transformed Cells Protein; Cell Proliferation-Inducing Gene 2 Protein; Hereditary Mixed Polyposis Syndrome; Colorectal Adenoma And Carcinoma 1; DAN Domain Family Member 2; Gremlin-1; CRAC1; IHG-2; Gremlin 1, Cysteine Knot Superfamily, Homolog (Xenopus Laevis); Gremlin 1, Cysteine Knot Superfamily, Homolog; Down-Regulated In Mos-Transformed Cells; Increased In High Glucose Protein; Increased In High Glucose-2; Gremlin 1-Like Protein; Gremlin 1; C15DUPq; DUP15q; CRCS4; HMPS1; MPSH; PIG2

Background

Gene Name: GREM1

NCBI Gene Entry: [26585](#)

UniProt Entry: [O60565](#)

Application Information

Molecular Weight: Predicted, 21 kDa

Clonality: Mouse polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC)

Immunogen

A synthesized peptide derived from human GREM1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Immunohistochemistry (IHC): 1:50-1:200

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Note: This product is for research use only.

Validation Data

