

Catalog #: 52927

Aliases

8-hydroxy-2'-deoxyguanosine; 2'-Deoxy-8-oxoguanosine; 7, 8-Dihydro-8-oxo-2'-deoxyguanosine; 8-Oxo-7, 8-dihydro-2'-deoxyguanosine; 8-Oxo-7, 8-dihydrodeoxyguanosine; 2'-Deoxy-8-oxo-D-guanosine; 8-Oxo-7, 8-dihydro-2 μ -deoxyguanosine, 8-Oxo-dG; 2'-Deoxy-8-hydroxyguanosine

Background

NCBI Gene Entry: [88847-89-6 \(CSA\)](#)

Application Information

Molecular Weight: Predicted, 0.283 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Species independent

Applications Tested: Immunohistochemistry (IHC), immunocytochemistry (IC)

Immunogen

8-OHdG

Isotype

Rabbit IgG

Storage Buffer

Supplied in 0.01M TBS (pH7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin300.

Recommended Dilutions

Immunohistochemistry (IHC): 1:100-1:500

Immunocytochemistry (IC): 1:100-1:500

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

