

Catalog #: 52929

Aliases

NOXO1; NADPH Oxidase Organizer 1; SH3PXD5; P41NOXA; P41NOXB; P41NOXC; SNX28; SH3 And PX Domain-Containing Protein 5; NADPH Oxidase Regulatory Protein; Nox-Organizing Protein 1; Nox Organizer 1; P41NOX; Regulatory Protein P41NOX

Background

Gene Name: NOXO1

NCBI Gene Entry: [124056](#)

UniProt Entry: [Q8NFA2](#)

Application Information

Molecular Weight: Predicted, 41 kDa; observed, 41 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NOXO1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of lysates from K562 cells, primary antibody was diluted at 1:1,000, 4°over night.