

Catalog #: 52930

Aliases

PERM1; PPARGC1 And ESRR Induced Regulator, Muscle 1; C1orf170; Peroxisome Proliferator-Activated Receptor Gamma Coactivator 1 And Estrogen-Related Receptor-Induced Regulator In Muscle 1; PGC-1 And ERR-Induced Regulator In Muscle Protein 1; PPARGC1 And ESRR-Induced Regulator In Muscle 1; RP11-54O7.8; MGC13275; PGC-1- And ERR-Induced Regulator In Muscle 1; PGC-1 And ERR-Induced Regulator In Muscle 1; Chromosome 1 Open Reading Frame 170

Background

Gene Name: PERM1

NCBI Gene Entry: [84808](#)

UniProt Entry: [Q5SV97](#)

Application Information

Molecular Weight: Predicted, 81 kDa; observed, 77 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human PERM1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,000

Note: This product is for research use only.

Validation Data

Western blotting analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1,000, 4°over night.