

## Catalog #: 5339

### Aliases

PFKP; Phosphofructokinase, Platelet; PFK-C; PFKF; ATP-Dependent 6-Phosphofructokinase, Platelet Type; Phosphofructo-1-Kinase Isozyme C; 6-Phosphofructokinase Type C; Phosphohexokinase; EC 2.7.1.11; ATP-PFK; PFK-P; 6-Phosphofructokinase, Platelet Type; Phosphofructokinase, Platelet Type; Phosphofructokinase 1; EC 2.7.1

### Background

Gene Name: PFKP

NCBI Gene Entry: [5214](#)

UniProt Entry: [Q01813](#)

### Application Information

Molecular Weight: Predicted, 86 kDa; observed, 75 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB1035

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

### Immunogen

Recombinant protein of human PFKP

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

---

#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

**Note:** This product is for research use only.

## Validation Data



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of PFKP expression in HepG2 cells using anti-PFKP antibody (Cat#5339, 1:2,000). Green, isotype control; red, PFKP.


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Immunocytochemical staining of HepG2 cells with anti-PFKP antibody (Cat#5339, 1:1,000). Nuclei were stained blue with DAPI; PFKP was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and Smart Gain: Medium. Scale bar, 20  $\mu$ m.


## SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

## ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-PFKP antibody (Cat#5339). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-PFKP antibody (Cat#5339, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**