

Catalog #: 5369

Aliases

AXL; AXL Receptor Tyrosine Kinase; UFO; JTK11; Tyro7; ARK; Tyrosine-Protein Kinase Receptor UFO; AXL Oncogene; EC 2.7.10.1; AXL Transforming Sequence/Gene; EC 2.7.10; AXL3

Background

Gene Name: AXL

NCBI Gene Entry: [558](#)

UniProt Entry: [P30530](#)

Application Information

Molecular Weight: Predicted, 98 kDa; observed, 138 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB1175

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human Axl

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-AXL antibody (Cat#5369). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-AXL antibody (Cat#5369, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM