

Catalog #: 5383

Aliases

RAB27A; RAB27A, Member RAS Oncogene Family; RAB27; HsT18676; RAM; GS2; Ras-Related Protein Rab-27A; GTP-Binding Protein Ram; Rab-27; Mutant Ras-Related Protein Rab-27A; EC 3.6.5.2

Background

Gene Name: RAB27A

NCBI Gene Entry: [5873](#)

UniProt Entry: [P51159](#)

Application Information

Molecular Weight: Predicted, 25 kDa; observed, 27 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB1245

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human RAB27

Isotype

Mouse IgG1 kappa

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western blotting (WB): 1:400-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-RAB27A antibody (Cat#5383). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-RAB27A antibody (Cat#5383, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).