

Catalog #: 5411

Aliases

DPP9; Dipeptidyl Peptidase 9; Dipeptidyl Peptidase-Like Protein 9; Dipeptidyl Peptidase IX; DPP IX; DPRP-2; DPLP9; DPRP2; DP9; Dipeptidyl Peptidase IV-Related Protein-2; Dipeptidyl Peptidase IV-Related Protein 2; Dipeptidylpeptidase 9; EC 3.4.14.5; IMD111; HATIS

Background

Gene Name: DPP9

NCBI Gene Entry: [91039](#)

UniProt Entry: [Q86TI2](#)

Application Information

Molecular Weight: Predicted, 98 kDa; observed, 98 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB1385

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

Recombinant protein of human DPP9

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Flow Cytometry (FCM): 1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Note: This product is for research use only.


Validation Data



Flow cytometric analysis of Dipeptidyl peptidase 9 expression in HepG2 cells using anti-Dipeptidyl peptidase 9 antibody (Cat#5411, 1:1,000). Green, isotype control; red, Dipeptidyl peptidase 9.



Immunocytochemical staining of HepG2 cells with anti-Dipeptidyl peptidase 9 antibody (Cat#5411, 1:1,000). Nuclei were stained blue with DAPI; Dipeptidyl peptidase 9 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.



Western blotting analysis using anti-dipeptidyl peptidase 9 antibody (Cat#5411). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-dipeptidyl peptidase 9 antibody (Cat#5411, 1:2,500) and HRP-conjugated goat anti-

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM