

Biotin-Conjugated Mouse Anti-Human IgG Fc Secondary Antibody



Catalog #: 543

Application Information

Clonality: Mouse Monoclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),
Immunocytochemistry (IC)

Conjugation: Biotin

Immunogen

Human IgG Fc

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.05% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Immunocytochemistry (IC): 1:500-1:1,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

