

# **Biotin-Conjugated Mouse Anti-Human IgG4 Secondary Antibody**



**Catalog #: 549**

## **Application Information**

Clonality: Mouse Monoclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),  
Immunocytochemistry (IC)

Conjugation: Biotin

## **Immunogen**

Human IgG4

## **Storage**

Store at -20 °C for one year.

## **Recommended Dilutions**

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

## **Validation Data**

---

### **SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### **ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**

