

Catalog #: 5638

Aliases

LY9; Lymphocyte Antigen 9; SLAMF3; CD229; MLY9; Hly9; Signaling Lymphocytic Activation Molecule 3; T-Lymphocyte Surface Antigen Ly-9; Cell Surface Molecule Ly-9; SLAM Family Member 3; CD229 Antigen

Background

Gene Name: LY9

NCBI Gene Entry: [4063](#)

UniProt Entry: [Q9HBG7](#)

Application Information

Molecular Weight: Predicted, 72 kDa; observed, 60-130 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB4640

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human LY9

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-CD229 antibody (Cat#5638). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-CD229 antibody (Cat#5638, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM