

Catalog #: 5643

Aliases

GMEB1; Glucocorticoid Modulatory Element Binding Protein 1; P96PIF; PIF96; Glucocorticoid Modulatory Element-Binding Protein 1; Parvovirus Initiation Factor P96; DNA-Binding Protein P96PIF; PIF P96; GMEB-1

Background

Gene Name: GMEB1

NCBI Gene Entry: [10691](#)

UniProt Entry: [Q9Y692](#)

Application Information

Molecular Weight: Predicted, 63 kDa; observed, 80 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB4665

Species Reactivity: Human, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM)

Immunogen

Recombinant protein of human GMEB1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-GMEB1 antibody (Cat#5643). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-GMEB1 antibody (Cat#5643, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of GMEB1 expression in HAP-1 cells using anti-GMEB1 antibody (Cat#5643, 1:2,000). Green, isotype control; red, GMEB1.