

Catalog #: 5657

Aliases

SEMA7A; Semaphorin 7A (JohnMiltonHagen Blood Group); CD108; H-Sema-L; SEMAL; Sema Domain, Immunoglobulin Domain (Ig), And GPI Membrane Anchor, (Semaphorin) 7A (JMH Blood Group); Semaphorin 7A, GPI Membrane Anchor (John Milton Hagen Blood Group); John-Milton-Hagen Human Blood Group Ag; JMH Blood Group Antigen; Semaphorin-7A; Semaphorin-K1; Semaphorin-L; Sema K1; CDw108; Sema L; Sema Domain, Immunoglobulin Domain (Ig), And GPI Membrane Anchor, (Semaphorin) 7A; Sema Domain, Immunoglobulin Domain (Ig), And GPI Membrane Anchor, 7A; Semaphorin 7A (John Milton Hagen Blood Group); John Milton Hagen Blood Group H-Sema K1; John Milton Hagen Blood Group; CD108 Antigen; H-Sema K1; H-SEMA-K1; PFIC11; SEMAK1; JMH

Background

Gene Name: SEMA7A

NCBI Gene Entry: [8482](#)

UniProt Entry: [O75326](#)

Application Information

Molecular Weight: Predicted, 75 kDa; observed, 75 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB4735

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human SEMA7A

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data



Western blotting analysis using anti-SEMA7A antibody (Cat#5657). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-SEMA7A antibody (Cat#5657, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM