

Catalog #: 5780

Aliases

STX17; Syntaxin 17; Syntaxin-17; FLJ20651

Background

Gene Name: STX17

NCBI Gene Entry: [55014](#)

UniProt Entry: [P56962](#)

Application Information

Molecular Weight: Predicted, 33 kDa; observed, 34 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB5350

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human STX17

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:400-1:2,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-syntaxin 17 antibody (Cat#5780). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-syntaxin 17 antibody (Cat#5780, 1:2,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM