

## Catalog #: 5782

## Aliases

TRIM47; Tripartite Motif Containing 47; RNF100; GOA; Gene Overexpressed In Astrocytoma Protein; Tripartite Motif-Containing Protein 47; E3 Ubiquitin-Protein Ligase TRIM47; RING Finger Protein 100; Tripartite Motif-Containing 47; EC 2.3.2.27

## Background

Gene Name: TRIM47

NCBI Gene Entry: [91107](#)

UniProt Entry: [Q96LD4](#)

## Application Information

Molecular Weight: Predicted, 70 kDa; observed, 74 kDa

Clonality: Mouse oligoclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

## Immunogen

Recombinant protein of human TRIM47

## Isotype

Mouse IgG

## Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis using anti-TRIM47 antibody (Cat#5782). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TRIM47 antibody (Cat#5782, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ<sup>TM</sup> ECL Substrate Kit (Cat#226).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)