

Catalog #: 5786

Aliases

HNRNPF; Heterogeneous Nuclear Ribonucleoprotein F; HNRPF; Nucleolin-Like Protein
Mcs94-1; HnRNP F Protein; OK/SW-Cl.23; Mcs94-1; HnRNP F

Background

Gene Name: HNRNPF
NCBI Gene Entry: [3185](#)
UniProt Entry: [P52597](#)

Application Information

Molecular Weight: Predicted, 46 kDa; observed, 46 kDa
Clonality: Mouse monoclonal antibody
Clone ID: 25GB5380
Species Reactivity: Human, mouse, rat
Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human HNRNPF

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:2,000-1:10,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-heterogeneous nuclear ribonucleoprotein F antibody (Cat#5786). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-heterogeneous nuclear ribonucleoprotein F antibody (Cat#5786, 1:10,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM