

Catalog #: 5790

Aliases

FAF1; Fas Associated Factor 1; UBXD12; UBXN3A; HFAP1; CGI-03; HFAP1s; Fas (TNFRSF6) Associated Factor 1; UBX Domain-Containing Protein 12; UBX Domain-Containing Protein 3A; TNFRSF6-Associated Factor 1; FAS-Associated Factor 1; UBX Domain Protein 3A

Background

Gene Name: FAF1

NCBI Gene Entry: [11124](#)

UniProt Entry: [Q9UNN5](#)

Application Information

Molecular Weight: Predicted, 74 kDa; observed, 74 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB5400

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human FAF1

Isotype

Mouse IgG1

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-fas associated factor 1 antibody (Cat#5790). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-fas associated factor 1 antibody (Cat#5790, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM