

## Catalog #: 5806

### Aliases

HNRNPH1; Heterogeneous Nuclear Ribonucleoprotein H1; Heterogeneous Nuclear Ribonucleoprotein H; HnRNPH; HNRPH1; Heterogeneous Nuclear Ribonucleoprotein H1 (H); HNRPH; Epididymis Secretory Sperm Binding Protein; HnRNP H; NEDCDS

### Background

Gene Name: HNRNPH1

NCBI Gene Entry: [3187](#)

UniProt Entry: [P31943](#)

### Application Information

Molecular Weight: Predicted, 49 kDa; observed, 49 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB5485

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

### Immunogen

Recombinant protein of human HNRNPH1

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

**Note:** This product is for research use only.

### Validation Data

---


#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis using anti-HNRNPH1 antibody (Cat#5806). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HNRNPH1 antibody (Cat#5806, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)