

# **FITC-Conjugated Mouse Anti-Human IgM Secondary Antibody**



**Catalog #: 581**

## **Application Information**

Clonality: Mouse Monoclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Immunohistochemistry (IHC), Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: FITC

## **Immunogen**

Human IgM

## **Storage Buffer**

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.05% sodium azide.

## **Storage**

Store at -20 °C for one year.

## **Recommended Dilutions**

Immunohistochemistry (IHC): 1:100-1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:1,000-1:4,000

## **Validation Data**

---

### **SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### **ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**

**FITC-Conjugated Mouse Anti-Human IgM  
Secondary Antibody**

PAGE 2



---

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**