

Catalog #: 5820

Aliases

GSTK1; Glutathione S-Transferase Kappa 1; GST13; GST Class-Kappa; GST 13-13; GSTK1-1; HGSTK1; Glutathione S-Transferase Subunit 13 Homolog; Glutathione S-Transferase Subunit 13; Glutathione S-Transferase K1; EC 2.5.1.18; GST13-13; GST;

Background

Gene Name: GSTK1

NCBI Gene Entry: [373156](#)

UniProt Entry: [Q9Y2Q3](#)

Application Information

Molecular Weight: Predicted, 25 kDa; observed, 25 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB5555

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Glutathione S Transferase kappa 1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis using anti-GSTK1 antibody (Cat#5820). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-GSTK1 antibody (Cat#5820, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQTM ECL Substrate Kit (Cat#716).