

Catalog #: 5823

Aliases

PDLIM7; PDZ And LIM Domain 7; ENIGMA; PDZ And LIM Domain 7 (Enigma); PDZ And LIM Domain Protein 7; Protein Enigma; LMP; Lim Mineralization Protein 3; LIM Mineralization Protein; LIM Domain Protein; 1110003B01Rik; LMP1; LMP3

Background

Gene Name: PDLIM7

NCBI Gene Entry: [9260](#)

UniProt Entry: [Q9NR12](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 49 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB5570

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

Recombinant protein of human PDLIM7

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-PDLIM7 antibody (Cat#5823). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-PDLIM7 antibody (Cat#5823, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM