

## Catalog #: 6011

### Aliases

VBP1; VHL Binding Protein 1; PFDN3; PFD3; Von Hippel-Lindau Binding Protein 1; Prefoldin Subunit 3; HIBBJ46; VBP-1; Von Hippel-Lindau-Binding Protein 1; VHL-Binding Protein 1; Prefoldin 3

### Background

Gene Name: VBP1

NCBI Gene Entry: [7411](#)

UniProt Entry: [P61758](#)

### Application Information

Molecular Weight: Predicted, 23 kDa; observed, 23 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 25GB6375

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

### Immunogen

Recombinant protein of human VBP1

### Isotype

Mouse IgG1

### Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

### Storage

Store at -20 °C for one year.

### Recommended Dilutions

Western Blotting (WB): 1:500-1:2,500

**Note:** This product is for research use only.

### Validation Data

---


#### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

#### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis using anti-VBP1 antibody (Cat#6011). Total cell lysates (30  $\mu$ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-VBP1 antibody (Cat#6011, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)