

Anti-FNTA Recombinant Rabbit Monoclonal Antibody



Catalog #: 6098

Aliases

FNTA; Farnesyltransferase, CAAX Box, Subunit Alpha; Protein Farnesyltransferase/ Geranylgeranyltransferase Type-1 Subunit Alpha; PGGT1A; PTAR2; FPTA; Protein Prenyltransferase Alpha Subunit Repeat Containing 2; Ras Proteins Prenyltransferase Subunit Alpha; Farnesyltransferase, CAAX Box, Alpha; GGTase-I-Alpha; FTase-Alpha; Type I Protein Geranyl-Geranyltransferase Alpha Subunit; Type I Protein Geranyl-Geranyltransferase Subunit Alpha; Farnesyl-Protein Transferase Alpha-Subunit; CAAX Farnesyltransferase Subunit Alpha; EC 2.5.1.58; EC 2.5.1.59

Background

Gene Name: FNTA

NCBI Gene Entry: [2339](#)

UniProt Entry: [P49354](#)

Application Information

Molecular Weight: Predicted, 44 kDa; observed, 44 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 25GB6855

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human FNTA

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-FNTA Recombinant Rabbit Monoclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



Western blotting analysis using anti-FNTA antibody (Cat#6098). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-FNTA antibody (Cat#6098, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using NaQTM ECL Substrate Kit (Cat#716).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM