

Catalog #: 615

Application Information

Clonality: Mouse Monoclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC)

Immunogen

Human IgG3

Storage Buffer

Supplied in PBS (pH 7.2), and 0.05% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

ELISA (E): 1:5,000-1:20,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Validation Data

