

Anti-Alpha Skeletal Muscle Actin Recombinant Rabbit Monoclonal Antibody



Catalog #: 6171

Aliases

Actin Alpha 1, Skeletal Muscle; NEM3; ACTA; Actin, Alpha Skeletal Muscle; Nemaline Myopathy Type 3; Alpha-Actin-1; EC 3.6.4.-; CMYP2A; CMYP2B; CMYP2C; CFTD1; CFTDM; ASMA; CFTD; MPFD; NEM1; NEM2; SHPM

Background

Gene Name: ACTA1
NCBI Gene Entry: [58](#)
UniProt Entry: [P68133](#)

Application Information

Molecular Weight: Predicted, 42 kDa; observed, 42 kDa
Clonality: Rabbit monoclonal antibody
Clone ID: 23GB1315
Species Reactivity: Human, mouse, rat
Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Actin

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000
Flow Cytometry (FCM): 1:2,000
Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Alpha Skeletal Muscle Actin Recombinant Rabbit Monoclonal Antibody


PAGE 2

Validation Data



Copyright ©2024 Genuin Biotechnologies LLC

Flow cytometric analysis of Actin expression in HT-1080 cells using Actin antibody (Cat#6171, 1:2,000). Green, isotype control; red, Actin.



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Alpha Skeletal Muscle Actin Recombinant Rabbit Monoclonal Antibody

PAGE 3

Immunocytochemical staining of HT-1080 cells with Actin antibody (Cat#6171, 1:1,000). Nuclei were stained blue with DAPI; Actin was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.



Western blotting analysis using anti-Actin antibody (Cat#6171). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Actin antibody (Cat#6171, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM