

KD-Validated Anti-TAR DNA binding protein Recombinant Rabbit Monoclonal Antibody



Catalog #: 61846

Aliases

TARDBP; TAR DNA Binding Protein; TDP-43; ALS10; TAR DNA-Binding Protein 43; TAR DNA-Binding Protein-43; TDP43

Background

Gene Name: TARDBP

NCBI Gene Entry: [23435](#)

UniProt Entry: [Q13148](#)

Application Information

Molecular Weight: Predicted, 45 kDa, observed, 40 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB4730

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human TDP43

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM




Flow cytometric analysis of TAR DNA binding protein expression in C2C12 cells using TAR DNA binding protein antibody (Cat#61846, 1:2,000). Green, isotype control; red, TAR DNA binding protein.



Immunocytochemical staining of C2C12 cells with TAR DNA binding protein antibody (Cat#61846, 1:1,000). Nuclei were stained blue with DAPI; TAR DNA binding protein was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.


KD-Validated Anti-TAR DNA binding protein Recombinant Rabbit Monoclonal Antibody

PAGE 3



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-TAR DNA binding protein antibody (Cat#61846). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TAR DNA binding protein antibody (Cat#61846, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-TAR DNA binding protein antibody (Cat#61846). TAR DNA binding protein expression in wild type (WT) and TAR DNA binding protein shRNA knockdown (KD) HeLa cells with 30 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

KD-Validated Anti-TAR DNA binding protein Recombinant Rabbit Monoclonal Antibody

PAGE 4

incubated with anti-TAR DNA binding protein antibody (Cat#61846, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM