

KD-Validated Anti-DPP9 Mouse Monoclonal Antibody



Catalog #: 63848

Aliases

DPP9; Dipeptidyl Peptidase 9; Dipeptidyl Peptidase-Like Protein 9; Dipeptidyl Peptidase IX; DPP IX; DPRP-2; DPLP9; DPRP2; DP9; Dipeptidyl Peptidase IV-Related Protein-2; Dipeptidyl Peptidase IV-Related Protein 2; Dipeptidylpeptidase 9; EC 3.4.14.5; IMD111; HATIS

Background

Gene Name: DPP9

NCBI Gene Entry: [91039](#)

UniProt Entry: [Q86TI2](#)

Application Information

Molecular Weight: Predicted, 98 kDa; observed, 98 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB8170

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

Recombinant protein of human DPP9

Isotype

Mouse IgG2b

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

KD-Validated Anti-DPP9 Mouse Monoclonal Antibody

PAGE 2

Note: This product is for research use only.

Validation Data



Dipeptidyl peptidase 9-Alexa Fluor® 647

Flow cytometric analysis of Dipeptidyl peptidase 9 expression in HepG2 cells using anti-Dipeptidyl peptidase 9 antibody (Cat#63848, 1:1,000). Green, isotype control; red, Dipeptidyl peptidase 9.


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Immunocytochemical staining of HepG2 cells with anti-Dipeptidyl peptidase 9 antibody (Cat#63848, 1:500). Nuclei were stained blue with DAPI; Dipeptidyl peptidase 9 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.




Western blotting analysis using anti-dipeptidyl peptidase 9 antibody (Cat#63848). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-dipeptidyl peptidase 9 antibody (Cat#63848, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using

KD-Validated Anti-DPP9 Mouse Monoclonal Antibody

PAGE 4

NaQ™ ECL Substrate Kit (Cat#716).



Western blotting analysis using anti-dipeptidyl peptidase 9 antibody (Cat#63848). Dipeptidyl peptidase 9 expression in wild type (WT) and dipeptidyl peptidase 9 (DPP9) shRNA knockdown (KD) HeLa cells with 20 µg of total cell lysates. β-Tubulin serves as a loading control. The blot was incubated with anti-dipeptidyl peptidase 9 antibody (Cat#63848, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM