

# KD-Validated Anti-TRIM47 Mouse Monoclonal Antibody



**Catalog #: 65138**

## Aliases

TRIM47; Tripartite Motif Containing 47; RNF100; GOA; Gene Overexpressed In Astrocytoma Protein; Tripartite Motif-Containing Protein 47; E3 Ubiquitin-Protein Ligase TRIM47; RING Finger Protein 100; Tripartite Motif-Containing 47; EC 2.3.2.27

## Background

Gene Name: TRIM47

NCBI Gene Entry: [91107](#)

UniProt Entry: [Q96LD4](#)

## Application Information

Molecular Weight: Predicted, 70 kDa; observed, 74 kDa

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15905

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB), flow cytometry (FCM)

## Immunogen

Recombinant protein of human TRIM47

## Isotype

Mouse IgG1

## Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:200-1:1,000

Flow Cytometry (FCM): 1:2,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS


SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)




Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of TRIM47 expression in H9c2 cells using anti-TRIM47 antibody (Cat#65138, 1:2,000). Green, isotype control; red, TRIM47.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-TRIM47 antibody (Cat#65138). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-TRIM47 antibody (Cat#65138, 1:1,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using NaQ™ ECL Substrate Kit (Cat#716).



Western blotting analysis using anti-TRIM47 antibody (Cat#65138). TRIM47 expression in wild-type (WT) and TRIM47 knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-TRIM47 antibody (Cat#65138, 1:1,000) and HRP-conjugated goat anti-mouse secondary antibody (Cat#101, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).