

KD-Validated Anti-Annexin A7 Recombinant Rabbit Monoclonal Antibody



Catalog #: 69581

Aliases

ANXA7; Annexin A7; ANX7; Annexin VII; Annexin-7 ; SNX; SYNEXIN; Synexin

Background

Gene Name: ANXA7

NCBI Gene Entry: [310](#)

UniProt Entry: [P20073](#)

Application Information

Molecular Weight: Predicted, 53 kDa, observed, 50 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 23GB1235

Species Reactivity: Human

Applications Tested: Western blotting (WB), flow cytometry (FCM), immunocytochemistry (IC)

Immunogen

A synthesized peptide derived from human Annexin A7

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Flow Cytometry (FCM): 1:2,000

Immunocytochemistry (IC): 1:100-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Flow cytometric analysis of Annexin A7 expression in HepG2 cells using Annexin A7 antibody (Cat#69581, 1:2,000). Green, isotype control; red, Annexin A7.



Immunocytochemical staining of HepG2 cells with Annexin A7 antibody (Cat#69581, 1:1,000). Nuclei were stained blue with DAPI; Annexin A7 was stained magenta with Alexa Fluor® 647.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

KD-Validated Anti-Annexin A7 Recombinant Rabbit Monoclonal Antibody

PAGE 3

Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar: 20 μ m.



Western blotting analysis using anti-Annexin A7 antibody (Cat#69581). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-Annexin A7 antibody (Cat#69581, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000) respectively. Image was developed using FeQ™ ECL Substrate Kit (Cat#226).

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

