

# AMCA-Conjugated Rabbit Anti-Mouse IgG (H&L) Secondary Antibody



**Catalog #: 750**

## Application Information

Clonality: Rabbit Polyclonal antibody

Species Reactivity: Mouse

Applications Tested: ELISA (E), Immunohistochemistry (IHC), Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: AMCA

## Immunogen

Mouse IgG

## Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Immunohistochemistry (IHC): 1:100-1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:1,000-1:4,000

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

**AMCA-Conjugated Rabbit Anti-Mouse IgG  
(H&L) Secondary Antibody**

PAGE 2



---

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**