

# AP-Conjugated Rabbit Anti-Rat IgG (H&L) Secondary Antibody



**Catalog #: 760**

## Application Information

Clonality: Rabbit Polyclonal antibody

Species Reactivity: Rat

Applications Tested: ELISA (E), Western Blotting (WB), Immunohistochemistry (IHC),  
Immunocytochemistry (IC)

Conjugation: AP

## Immunogen

Rat IgG

## Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

ELISA (E): 1:10,000-1:50,000

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:500-1:1,000

Immunocytochemistry (IC): 1:500-1:1,000

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

**AP-Conjugated Rabbit Anti-Rat IgG (H&L)  
Secondary Antibody**

PAGE 2



---

**SUPPORT**

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

**ORDERS**

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

**WWW.GENUINBIOTECH.COM**