

Texas Red-Conjugated Goat Anti-Human IgG Fab Secondary Antibody



Catalog #: 846

Application Information

Clonality: Rabbit Polyclonal antibody

Species Reactivity: Human

Applications Tested: ELISA (E), Immunohistochemistry (IHC), Immunocytochemistry (IC), Flow Cytometry (FCM)

Conjugation: Texas Red

Immunogen

Human IgG Fab

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, 1% BSA and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Immunohistochemistry (IHC): 1:100-1:1,000

Immunocytochemistry (IC): 1:100-1:1,000

Flow Cytometry (FCM): 1:1,000-1:4,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

**Texas Red-Conjugated Goat Anti-Human IgG
Fab Secondary Antibody**

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM