

HRP-Conjugated Mouse Anti Rabbit IgG (Light chain specific) Secondary Antibody



Catalog #: 881

Application Information

Clonality: Mouse Monoclonal antibody

Species Reactivity: Rabbit

Applications Tested: Western Blotting (WB)

Conjugation: HRP

Immunogen

Rabbit IgG light chain

Isotype

Mouse IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:10,000-1:50,000

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

