

Rabbit Anti Sheep/Goat IgG (H&L) Secondary Antibody



Catalog #: 894

Application Information

Clonality: Rabbit polyclonal antibody

Species Reactivity: Sheep/Goat

Applications Tested: Western blotting (WB), immunohistochemistry (IHC), ELISA (E)

Immunogen

Sheep/Goat IgG

Storage Buffer

Supplied in PBS (pH 7.2) containing 50% glycerol, and 0.05% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:5,000-1:20,000

Immunohistochemistry (IHC): 1:200-1:1,000

ELISA (E): 1:5,000-1:20,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

