

# Human CDKN1B Knockdown Cell Line (WB-Validated)



**Catalog #: C1241**

## Aliases

CDKN1B; Cyclin Dependent Kinase Inhibitor 1B; KIP1; P27KIP1; Cyclin-Dependent Kinase Inhibitor 1B (P27, Kip1); Cyclin-Dependent Kinase Inhibitor 1B; Cyclin-Dependent Kinase Inhibitor P27; P27Kip1; CDKN4; MEN1B; MEN4; P27

## Background

Gene Name: CDKN1B  
NCBI Gene Entry: [1027](#)

## Storage

Store at liquid nitrogen for 1 year.

## Kit Components

1. Human CDKN1B Knockdown Cell Line (Wb-Validated)
2. Wild-type cell line

## Parental Cell Line

Human cell line supplied by the client

## Validation Methods

RT-qPCR, Western blotting (WB)

## Shipping

Shipped on Dry Ice.

## Instructions For Use

This knockdown cell line should be paired with wild-type cell line for use.

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT


SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041



[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

Human CDKN1B Knockdown Cell Line (WB-Validated)



| Genotype                      | Ct Value |
|-------------------------------|----------|
| Wild-Type                     | 27.57    |
| Knock-Down                    | 29.16    |
| $\Delta Ct (Ct_{KD}-Ct_{WT})$ | 1.59     |
| % mRNA Reduction              | ↓ 67%    |

RT-qPCR analysis. HeLa cells were infected with CDKN1B-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers.  $\Delta Ct (Ct_{KD}-Ct_{WT})$  was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula:  $(1-1/2^{\Delta Ct}) \times 100\%$ .



SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

Western blotting analysis. CDKN1B protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting.  $\beta$ -Tubulin served as a loading control. The blots were incubated with primary antibodies (Cat#69245, 1:5,000) against CDKN1B and  $\beta$ -Tubulin, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000). Images were developed using FeQ™ ECL Substrate Kit (Cat#226).