

# Human SOD2 Knockdown Cell Line (WB-Validated)



**Catalog #: C62706**

## Aliases

SOD2; Superoxide Dismutase 2; GCInc1; MnSOD; IPOB; GC1; Superoxide Dismutase [Mn], Mitochondrial; Superoxide Dismutase 2, Mitochondrial; EC 1.15.1.1; Epididymis Secretory Sperm Binding Protein; Manganese-Containing Superoxide Dismutase; Gastric Cancer–Associated LncRNA 1; Gastric Cancer-Associated LncRNA 1; Manganese Superoxide Dismutase; Mangano-Superoxide Dismutase; Mn Superoxide Dismutase; Indophenoloxidase B; Mn-SOD; IPO-B; MVCD6

## Background

Gene Name: SOD2

NCBI Gene Entry: [6648](#)

## Storage

Store at liquid nitrogen for 1 year.

## Kit Components

1. Human SOD2 Knockdown Cell Line (Wb-Validated)
2. Wild-type cell line

## Parental Cell Line

Human cell line supplied by the client

## Validation Methods

RT-qPCR, Western blotting (WB)

## Shipping

Shipped on Dry Ice.

## Instructions For Use

This knockdown cell line should be paired with wild-type cell line for use.

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)


## Human SOD2 Knockdown Cell Line (WB-Validated)

PAGE 2



Genotype	Ct Value
Wild-Type	16.07
Knock-Down	21.50
$\Delta Ct (Ct_{KD} - Ct_{WT})$	5.43
% mRNA Reduction	↓ 98%

RT-qPCR analysis. HeLa cells were infected with SOD2-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers.  $\Delta Ct (Ct_{KD} - Ct_{WT})$  was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula:  $(1 - 1/2^{\Delta Ct}) \times 100\%$ .



Western blotting analysis. SOD2 protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting.  $\beta$ -Tubulin served as a loading control. The blots were incubated with primary antibodies against SOD2 and  $\beta$ -Tubulin, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM