

Human GLS Knockdown Cell Line (WB-Validated)



Catalog #: C63470

Aliases

GLS; Glutaminase; GLS1; K-Glutaminase; KIAA0838; GAC; GAM; KGA; Glutaminase Kidney Isoform, Mitochondrial; L-Glutamine Amidohydrolase; Glutaminase C; EC 3.5.1.2; Glutaminase, Phosphate; Activated; CASGID; EIEE71; AAD20; DEE71; GDPAG

Background

Gene Name: GLS

NCBI Gene Entry: [2744](#)

Storage

Store at liquid nitrogen for 1 year.

Kit Components

1. Human GLS Knockdown Cell Line (Wb-Validated)
2. Wild-type cell line

Parental Cell Line

Human cell line supplied by the client

Validation Methods

RT-qPCR, Western blotting (WB)

Shipping

Shipped on Dry Ice.

Instructions For Use

This knockdown cell line should be paired with wild-type cell line for use.

Note: This product is for research use only.

Validation Data

SUPPORT


SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM


Human GLS Knockdown Cell Line (WB-Validated)



Genotype	Ct Value
Wild-Type	20.35
Knock-Down	30.3
$\Delta Ct (Ct_{KD}-Ct_{WT})$	9.95
% mRNA Reduction	↓ 99.9%

Copyright ©2024 Genuin Biotechnologies LLC

RT-qPCR analysis. HeLa cells were infected with GLS-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers. $\Delta Ct (Ct_{KD}-Ct_{WT})$ was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula: $(1-1/2^{\Delta Ct}) \times 100\%$.



Copyright ©2024 Genuin Biotechnologies LLC

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Western blotting analysis. GLS protein expression in wild-type (WT) and shRNA knockdown (KD) HeLa cells was detected using Western blotting. β -Tubulin served as a loading control. The blots were incubated with primary antibodies (Cat#63470, 1:5,000) against GLS and β -Tubulin, respectively, followed by incubating with HRP-conjugated goat anti-rabbit secondary antibody (Cat#201, 1:20,000). Images were developed using FeQ™ ECL Substrate Kit (Cat#226).