

Human RBX1 Knockdown Cell Line (WB-Validated)



Catalog #: C63745

Aliases

RBX1; Ring-Box 1; RNF75; ROC1; Regulator Of Cullins 1; BA554C12.1; E3 Ubiquitin-Protein Transferase RBX1; E3 Ubiquitin-Protein Ligase RBX1; RING Finger Protein 75; RING-Box Protein 1; Ring-Box 1, E3 Ubiquitin Protein Ligase; ZYP Protein; Protein ZYP; EC 2.3.2.27; EC 2.3.2.32; Rbx1

Background

Gene Name: RBX1

NCBI Gene Entry: [9978](#)

Storage

Store at liquid nitrogen for 1 year.

Kit Components

1. Human RBX1 Knockdown Cell Line (Wb-Validated)
2. Wild-type cell line

Parental Cell Line

Human cell line supplied by the client

Validation Methods

RT-qPCR, Western blotting (WB)

Shipping

Shipped on Dry Ice.

Instructions For Use

This knockdown cell line should be paired with wild-type cell line for use.

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



RT-qPCR analysis. HT-1080 cells were infected with RBX1-specific shRNA lentiviral particles, total RNA was extracted from wild-type and knockdown cells, RT-qPCR was performed using gene-specific primers. $\Delta Ct (Ct_{KD} - Ct_{WT})$ was used to calculate mRNA reduction (%) between wild-type and knockdown cells using the following formula: $(1 - 1/2^{\Delta Ct}) \times 100\%$.



Western blotting analysis. RBX1 protein expression in wild-type (WT) and shRNA knockdown (KD) HT-1080 cells was detected using Western blotting. Hsp90 α served as a loading control. The blots were incubated with primary antibodies against RBX1 and Hsp90 α , respectively, followed by incubating with HRP-conjugated goat anti-mouse secondary antibody. Images were developed using FeQ™ ECL Substrate Kit.